J Antimicrob Chemother 2023; **78**: 2770–2771 https://doi.org/10.1093/jac/dkad243 Advance Access publication 7 August 2023

Rapid detection of temocillin resistance in Enterobacterales

Jacqueline Findlay¹*, Laurent Poirel (1)^{1,2,3} and Patrice Nordmann^{1,2,3,4}

¹Medical and Molecular Microbiology, Faculty of Science and Medicine, University of Fribourg, Chemin du Musée 18, CH-1700 Fribourg, Switzerland; ²Swiss National Reference Center for Emerging Antibiotic Resistance (NARA), University of Fribourg, Fribourg, Switzerland; ³INSERM European Unit (IAME, France), University of Fribourg, Fribourg, Switzerland; ⁴Institute for Microbiology, University of Lausanne and University Hospital Centre, Lausanne, Switzerland

*Corresponding author. E-mail: jacqueline.findlay@unifr.ch

Temocillin is a semi-synthetic $6-\alpha$ -methoxy derivative of ticarcillin, first developed in 1981.¹ Despite being developed over 40 years ago, temocillin's use as an antimicrobial agent was largely overlooked as a treatment option for infections caused by Gram-negative bacteria due to its poor activity against nonfermenters, including *Pseudomonas* spp. and *Acinetobacter* spp.¹ Temocillin has been demonstrated to have an affinity for penicillinbinding protein 3 in Escherichia coli² and remarkable stability against a plethora of β-lactamases including ESBLs and AmpCs (both plasmid and chromosomal);^{3,4} however, it is not clinically useful against bacteria producing class B (e.g. NDM) or class D (e.g. OXA-48-type) carbapenemases, because those latter enzymes readily hydrolyse this antibiotic, often leading to high MICs of temocillin.⁵ A recent study evaluating the use of temocillin against ESBL- and AmpC-producing Enterobacterales confirmed the excellent activity of temocillin, and the few high MICs observed were linked mostly to the carriage of multiple β -lactamases by corresponding isolates rather than any single β -lactamase type.⁶ In recent years, a number of countries in Europe, including Belgium, France, Luxembourg and the UK, have revived the use of this antibiotic, predominantly for the treatment of urinary tract infections but also for bloodstream and lower respiratory tract infections.^{7,8} Subsequently, the surveillance of resistance to temocillin is essential and can be easily determined using either routine broth microdilution or disc diffusion testing; however, such tests require 18-24 h to achieve a result. In this study, we designed and evaluated a rapid (\sim 4 h) test that can determine resistance or susceptibility to temocillin in Enterobacterales.

A total of 100 clinical Enterobacterales isolates were tested, comprising 55 E. coli, 33 Klebsiella pneumoniae, 6 Enterobacter

cloacae, 4 Klebsiella aerogenes and 2 Klebsiella oxytoca isolates (Table S1, available as Supplementary data at JAC Online). Isolates were submitted to the Swiss National Reference Center for Emerging Antibiotic Resistance (NARA) from hospitals and clinics throughout Switzerland, and speciated by either MALDI-TOF MS (Bruker Microflex LT, Bruker Daltonik GmbH, Bremen, Germany) or API-20E tests (bioMérieux, https://www. biomerieux.com). Susceptibility testing was performed by broth microdilution according to CLSI guidelines,⁹ and resistance was determined according to EUCAST breakpoints (R > 16 mg/L),¹⁰ because there are no available guidelines for temocillin resistance according to CLSI. This identified that 38 tested isolates were temocillin susceptible, with MICs ranging from 2 to 16 mg/L, and the remaining 62 were resistant, with MICs ranging from 32 to >256 mg/L. The β -lactamase content of the tested isolates was determined by multiplex PCRs targeting CTX-M, CMY, DHA, SHV, TEM and OXA-1-like encoding genes,¹¹ and the carbapenemase gene content had previously been determined by PCR, and subseauent Sanaer sequencing of amplicons by NARA.¹²

The Rapid Temocillin NP test solution was prepared by mixing 6.25 g of cation-adjusted Mueller-Hinton (MH) broth (CA-MHB) powder (Bio-Rad), 0.0125 g of phenol red (Sigma Aldrich) and 225 mL of distilled water. The pH was adjusted to 7.3 by adding drops of 1 mol/L hydrochloric acid prior to autoclaving. After cooling to room temperature, 25 mL of a 10% anhydrous D(+)-alucose (Roth, Karlsruhe, Germany) (sterilized by filtration) was added, resulting in the final concentrations of 2.5% CA-MHB, 0.005% phenol red indicator and 1% glucose. The temocillin test solution was prepared with 21.3 mg/L of temocillin and 150 µL was added into one well of a 96-well polystyrene plate; 150 µL of test solution without temocillin was added into a second well. Inoculums were prepared from overnight cultures, grown on URISelect 4 agar (Bio-Rad, https://www.bio-rad.com), as follows: a 0.5 MacFarland suspension was prepared then diluted 1/10 in fresh MH broth; 50 µL of this suspension was then used to inoculate the plates (approx. 1.5×10^7 bacterial cells per well), resulting in a final temocillin concentration of 16 mg/L, before incubation at 37°C for 4 h. Negative control wells, both with and without temocillin, were inoculated with 50 µL 0.85% NaCl, diluted 1/10 in fresh MH broth, instead of the bacterial suspension. After 3 h and 4 h, plates were read and a colour change from red to yellow was interpreted as positive (resistant), confirming acid production via the metabolism of glucose in the test solution, whereas red was interpreted as negative (susceptible) (Figure 1).

The test presented a single false-negative result, this being an *E. coli* isolate (T32) harbouring bla_{TEM} as a sole acquired ß-lactamase gene, exhibiting an MIC of temocillin at 32 mg/L, therefore corresponding to a 100% specificity and 98.4% sensitivity. Most of the isolates exhibiting high-level temocillin

[©] The Author(s) 2023. Published by Oxford University Press on behalf of British Society for Antimicrobial Chemotherapy. This is an Open Access article distributed under the terms of the Creative Commons Attribution License (https://creativecommons.org/licenses/ by/4.0/), which permits unrestricted reuse, distribution, and reproduction in any medium, provided the original work is properly cited. 2770




Figure 1. An example of a negative and positive test result. 0 = test solution containing no temocillin; 16 = test solution containing temocillin at a final concentration of 16 mg/L. T17, *E. coli* isolate with an MIC of 16 mg/L, T72, *K. pneumoniae* isolate with an MIC of 64 mg/L. This figure appears in colour in the online version of *JAC* and in black and white in the print version of *JAC*.

resistance (MIC \geq 64 mg/L) resulted in a positive test result after only 3 h; however, those with resistant MICs closer to the breakpoint (32 mg/L) required 4 h. The test was shown to be effective regardless of the Enterobacterales species, β -lactamase content and overall mechanism of resistance. This rapid test could be easily implemented in a clinical laboratory and can be set up alongside routine antimicrobial susceptibility testing methodologies (AST), but providing a result 14–20 h earlier than traditional AST and potentially sparing the use of other β -lactams such as the carbapenems. These results show that as the use of temocillin to treat Gram-negative infections becomes more commonplace, such a test can prove useful in determining targeted rather than simply empirical therapy in a relatively short time frame.

Funding

This work was financed by the Swiss National Science Foundation (grant FNS 310030_1888801).

Transparency declarations

None to declare.

Supplementary data

Table S1 is available as Supplementary data at JAC Online.

References

1 Jules K, Neu HC. Antibacterial activity and beta-lactamase stability of temocillin. *Antimicrob Agents Chemother* 1982; **22**: 453–60. https://doi. org/10.1128/AAC.22.3.453

2 Labia R, Baron P, Masson JM *et al.* Affinity of temocillin for *Escherichia coli* K-12 penicillin-binding proteins. *Antimicrob Agents Chemother* 1984; **26**: 335–8. https://doi.org/10.1128/AAC.26.3.335

3 Livermore DM, Hope R, Fagan EJ *et al.* Activity of temocillin against prevalent ESBL- and AmpC-producing Enterobacteriaceae from southeast England. *J Antimicrob Chemother* 2006; **57**: 1012–4. https://doi. org/10.1093/jac/dkl043 **4** Kuch A, Zieniuk B, Żabicka D *et al.* Activity of temocillin against ESBL-, AmpC-, and/or KPC-producing Enterobacterales isolated in Poland. *Eur J Clin Microbiol Infect Dis* 2020; **39**: 1185–91. https://doi.org/10.1007/s10096-020-03844-5

5 Livermore DM, Warner M, Mushtaq S *et al.* What remains against carbapenem-resistant Enterobacteriaceae? Evaluation of chloramphenicol, ciprofloxacin, colistin, fosfomycin, minocycline, nitrofurantoin, temocillin and tigecycline. *Int J Antimicrob Agents* 2011; **37**: 415–9. https://doi. org/10.1016/j.ijantimicag.2011.01.012

6 Stewart AG, Henderson A, Bauer MJ *et al.* Activity of temocillin against third-generation cephalosporin-resistant *Escherichia coli* and *Klebsiella pneumoniae* bloodstream isolates from a clinical trial. *JAC Antimicrob Resist* 2021; **4**: dlab192. https://doi.org/10.1093/jacamr/dlab192

7 Livermore DM, Tulkens PM. Temocillin revived. *J Antimicrob Chemother* 2009; **63**: 243–5. https://doi.org/10.1093/jac/dkn511

8 Giske CG. Contemporary resistance trends and mechanisms for the old antibiotics colistin, temocillin, fosfomycin, mecillinam and nitrofurantoin. *Clin Microbiol Infect* 2015; **21**: 899–905. https://doi.org/10.1016/j.cmi. 2015.05.022

9 CLSI. Methods for Dilution Antimicrobial Susceptibility Tests for Bacteria That Grow Aerobically; Approved Standard—Tenth Edition: M07-A11. 2018.

10 European Committee on Antimicrobial Susceptibility Testing. Clinical breakpoints v13.0. https://www.eucast.org/fileadmin/src/media/PDFs/EUCAST_files/Breakpoint_tables/v_13.0_Breakpoint_Tables.pdf

11 Findlay J, Gould VC, North P et al. Characterization of cefotaxime-resistant urinary *Escherichia coli* from primary care in South-West England 2017–18. J Antimicrob Chemother 2020; **75**: 65–71. https://doi.org/10.1093/jac/dkz397

12 Poirel L, Walsh TR, Cuvillier V *et al.* Multiplex PCR for detection of acquired carbapenemase genes. *Diagn Microbiol Infect Dis* 2011; **70**: 119–23. https://doi.org/10.1016/j.diagmicrobio.2010.12.002

J Antimicrob Chemother 2023; **78**: 2771–2774 https://doi.org/10.1093/jac/dkad253 Advance Access publication 11 August 2023

Evaluation of the RESIST ACINETO multiplex immunochromatographic assay for detection of OXA-23-like, OXA-40/58-like and NDM carbapenemase production in *Acinetobacter baumannii*

Stefano Mancini¹*, Helena M. B. Seth-Smith¹, Natalia Kolesnik-Goldmann¹, Vladimira Hinic¹, Tim Roloff¹, Frank Imkamp¹ and Adrian Egli ¹

¹Institute of Medical Microbiology, University Zurich, Gloriastrasse 28/30, 8006 Zurich, Switzerland

*Corresponding author. E-mail: smancini@imm.uzh.ch

© The Author(s) 2023. Published by Oxford University Press on behalf of British Society for Antimicrobial Chemotherapy. This is an Open Access article distributed under the terms of the Creative Commons Attribution License (https://creativecommons.org/licenses/ by/4.0/), which permits unrestricted reuse, distribution, and reproduction in any medium, provided the original work is properly cited.